

Materials in PendoTECH® Single Use Pressure Sensors™


Introduction

The PendoTECH Single Use Pressure Sensors have proven to be an invaluable tool to measure pressure in many process operations. As more customers are in the process of qualifying them for use in their process or are considering using the sensors in certain processes, the information on the material used in the sensors becomes a required piece of information. This Technical Note is designed to provide some of the details on the materials.



Flow Path Materials

1. Sensor body (> 99% of surface area)-
Two materials are offered as alternatives for different customer applications, polycarbonate and polysulfone. Polysulfone offers better chemical compatibility including compatibility with NaOH (Application Note available).
 - a. Part numbers starting with PRESS
 - i. Covestro Makrolon Rx1805 polycarbonate
 - b. Part numbers starting with PREPS
 - i. Solvay UDEL P1700 polysulfone
2. Pressure Sensing Chip
 - a. The material that isolates the liquid in the fluid path and away from the electronics and serves to communicate pressure between the flow path and pressure measuring diaphragm is a two-component, cured in place insoluble dielectric silicone
 - i. Proprietary formulation
 - b. Molded plastic component that is affixed to the pressure sensor body, part of which is a shallow well to hold the insoluble silicone
 - i. Proprietary polycarbonate plastic formulation (tested for NaOH resistance by PendoTECH)
3. Adhesive
 - a. Used to affix the molded plastic component of the pressure sensing chip to the pressure sensor body
 - i. Dymax proprietary formulation



Component	Biocompatibility Testing ⁽¹⁾	Animal Derived Material Content ⁽²⁾
Covestro Makrolon Rx 1805	USP Class VI ISO 10993-5	Animal Derived Component Free
Solvay UDEL P1700	USP Class VI ISO 10993-5	Animal Derived Component Free
Dielectric Silicone	USP Class VI ISO 10993-5	Animal Derived Component Free
Pressure Sensing Chip Proprietary Plastic	USP Class VI ISO 10993-5	Animal Derived Component Free
Dymax Proprietary Adhesive	USP Class VI ISO 10993-5	Animal Derived Component Free

⁽¹⁾ Testing by PendoTECH

⁽²⁾ Letters on file from original manufacturer at PendoTECH