

WORKSHEET FOR:

PendoTECH Cable with RJ11 Phone Connector Adapter
For Single Use Pressure Sensor
Part# PDKT-650-XXX

RETURN VIA FAX +1-609-784-7889 or EMAIL: request@pendotech.com


ATTENTION: Custom Cables Representative

FROM: Name: _____
Company: _____
Phone: _____
Email: _____

PLEASE COMPLETE #1, #2, and #3 BELOW

#1 Electrical Specifications Existing Monitoring System/Sensor:
Input Impedance 270 Ohms or greater: YES or NO
Output Impedance 270 Ohms or greater: YES or NO
Excitation Voltage of 2.5 to 10 Volts DC: YES or NO

Clip Side



#2 Indicate desired length of cable (12 feet or shorter):
_____ feet

- Pin 4**
Signal +
Signal -
Excitation +
Excitation -
- Pin 3**
Signal +
Signal -
Excitation +
Excitation -
- Pin 2**
Signal +
Signal -
Excitation +
Excitation -
- Pin 1**
Signal +
Signal -
Excitation +
Excitation -

#3 Circle the function of the pin in diagram of your existing sensor (each selection can be made only once)

IMPORTANT NOTE: After completing #1 and #2 a sensor can be sent to PendoTECH for testing to confirm pin arrangement and to confirm the signal output (accuracy) is the same.