

PressureMAT® DPG

- Reads PendoTECH Single Use Pressure Sensors™
- Convenient handheld design
- Battery operated
- No calibration required when used with the pressure sensors (the standardized sensor output signal is read)
- One-touch button for sensor tare to atmospheric pressure
- Minimum / maximum pressure display buttons
- Comes with 12foot re-usable cable



PressureMAT DPG Specifications

- Reads sensors to their specification of 0-6psi +/- 2% of reading; 6-30psi +/- 3% of reading
- Uses 9V battery (approximate 14 hours of power)
- Auto shut-off after 10 minutes
- Resolution displays x.xx psi with updates >2X per second

Features of PendoTECH Pressure Sensors

- They come in a range of sizes including a luer fitting, hosebarb design, and sanitary flange fittings (as shown in picture to the right)
Click here for Pressure Sensor Data Sheet
- Unobstructed flow path and no dead legs
- Qualified for use with no failure to -11.5 to 75psi
- Compatible with gamma irradiation
- All materials meet USP Class VI
- Integral 12inch cable connects to reusable cable via molded connector
- Each sensor is tested for quality and accuracy during manufacturing and rejected if out of specification (NO CALIBRATION REQUIRED!)
- Manufactured in an ISO Class 7 clean room and is FDA registered, ISO 13485 facility
- Low cost
- Certificate of quality included with lot certification; individual NIST Calibration Certificates are optional
- May be non-invasively tested in-place via test port



Ordering Information

PMAT-DPG	Handheld monitor for one PendoTECH Pressure Sensor
----------	--

For warranty information see our website at <http://www.pendotech.com/warranty>
PendoTECH® is a registered trademark of PendoTECH, all rights reserved.