

PRODUCT PREVIEW - PLANNED AVAILABILITY 2022 Q3

Single Use In-line pH Sensor

Bringing proven and reliable METTLER TOLEDO single-use pH technology to an in-line format ideal for downstream processing applications



1/4 inch Hosebarb



3/4 inch Sanitary Flange

Key Features:

- Single-use pH sensor combining METTLER TOLEDO InSUS 307 pH probe technology with a single use flow cell designed by PendoTECH
- Designed for applications where in-line sensing is necessary
- Two sizes offered as a pre-assembled product
- Compatible with gamma irradiation
- No process calibration required
- Designed for use with existing METTLER TOLEDO transmitters

